

# Codonics Infinity 1.3.1 Operating Software

**Release Notes**

## **Version 1.3.1 Release Notes**

### **Summary**

Version 1.3.1 software corrects defects for the Codonics Infinity® C Medical Image Server. The software is compatible with all Infinity C and LX systems.

The Infinity Medical Image Server is an expandable DICOM storage appliance that provides a solution for storing and retrieving medical images. Applications include temporary storage of large DICOM data sets directly from any modality, and storage of digital images in an environment without a centralized, digital image archive.

### **Defects Corrected**

- ◆ **The Infinity C system clock resets to 1/1/2002 causing the Infinity WebViewer to not work.** When upgrading with 1.3.1 software, the software will correct Infinity C's time if it has already reset to 1/1/2002. The clock will be automatically reset to 6/1/2009. A user will then need to reset the clock and date to the proper time.

There is no effect to studies that are stored if the clock does change because the database doesn't store studies based on the system clock.

**NOTE:** During Infinity C system start-up with 1.3.0 or previous software versions, the red LEDs on the drives would remain lit. With 1.3.1 software and 105D BIOS, the red LEDs no longer remain lit during start-up.

**NOTE:** It is recommended to periodically check your date and time.

- ◆ **Software cannot be upgraded on Infinity C.** With 1.3.0 or previous software versions, a user may not be able to upgrade their software. Now, 1.3.1 software will allow a user to upgrade their software on Infinity C.

### **Technical Support**

If problems occur during software installation, contact Codonics Technical Support between the hours of 8:30AM and 5:30PM EST (weekends and U.S. holidays excluded).

Phone: +1-440-243-1198  
Email: support@codonics.com  
Website: www.codonics.com

*Get it all with just one call*  
**1-800-444-1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice.

Copyright © 2009 Codonics, Inc. Printed in the U.S.A. Part No. 901-223-001.03.



17991 Englewood Drive  
Cleveland, OH 44130 USA  
(440) 243-1198  
(440) 243-1334 Fax  
Email info@codonics.com  
www.codonics.com

Codonics Limited KK  
New Shibaura Bldg. F1  
1-3-11, Shibaura  
Minato-ku, Tokyo, 105-0023 JAPAN  
Phone: 81-3-5730-2297  
Fax: 81-3-5730-2295