

Infinity External USB Hard Drive Archiving


Overview

This technical brief explains how to archive studies from a Codonics® Infinity® Medical Image Server to an externally connected USB hard drive.

Archiving Studies

The following instructions outline the process for configuring an external USB hard drive and for exporting studies to it from an Infinity server.

1. Connect a USB hard drive, which is 2 TB or smaller, to the Infinity server.
2. Open the Study Manager for LX or C servers (or the **File Manager** for RDSS servers) and select **Options**, then **DICOM DIR Settings**.



3. Select the **Connections Configuration** tab in the **DICOM DIR Configuration** window and choose **Add**.
4. Select the **double arrow icon** to the right of the **DICOM DIR folder** entry field.



5. Using the **Browse for Folder** window, locate and select the USB hard drive and click the **OK** button.



NOTE: Optionally, select the **Scan Folder and Subfolders** check box to scan the files in addition to the DICOM DIR for later study searching on the USB HDD.

NOTE: The Windows drive letter assignment may change if any other USB device is connected to the Infinity server after configuration of the USB HDD.


6. Enter a **Configuration Name** for the USB hard drive on the **DICOM DIR Configuration** window.




7. If the USB hard drive has not had any studies previously sent to it from the Infinity, a prompt to enter a DICOM DIR name will be presented.
8. Enter a DICOM DIR name with a maximum of 16 characters.



9. The new USB hard drive will be presented as a location in the left window pane of Study/File Manager.
10. To send one or more studies to the USB hard drive, select each and choose **Export Study** by right-clicking.



11. On the **Connection Selection** window, choose the USB hard drive and then click the OK button.



USB Hard Drive Uses

The configured USB hard drive can now be used for the following user activities:

- ◆ Export of additional studies.
- ◆ Nearline storage access of the stored studies while connected to the Infinity.
- ◆ Offline storage access of the stored studies when connected to other Infinity servers or DICOM workstations.

Technical Support

Please contact Codonics Technical Support at any time.

Phone: +1.440.243.1198
 Email: support@codonics.com
 Website: www.codonics.com

Get it all with just one call
800.444.1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2014 Codonics, Inc. Printed in the U.S.A. Part No. 901-588-001.01.



17991 Englewood Drive
 Middleburg Heights, OH
 44130 USA
 +1.440.243.1198
 +1.440.243.1334 Fax
 Email info@codonics.com
 www.codonics.com

Codonics Trading Co, Ltd.
 317 Xianxia Rd. Building B
 Unit 1412
 Changning Dist., Shanghai
 P.R. China, 200051
 86-21-62787701
 86-21-62787719 Fax

Codonics Limited KK
 AQUACITY 9F,
 4-16-23, Shibaura
 Minato-ku, Tokyo,
 108-0023 JAPAN
 81-3-5730-2297
 81-3-5730-2295 Fax