

Codonics Virtua *syngo Media Viewer*

Technical Brief

Overview

The Codonics Virtua[®] Medical Disc Publisher provides support for Siemens *syngo*[®] Media Viewer as an optional DICOM Viewer, allowing the display of fused pre-aligned PET/CT or SPECT/CT images as well as stand-alone CT, MR, SPECT or PET studies. If the anatomical and functional data are acquired from different scanners, a workstation is required to align and export the fused data. *syngo* Media Viewer requires a Feature Key on Virtua to activate the Viewer for installation on the SmartDrive.

Minimum System Requirements

- ◆ Windows XP or Vista (Home and Professional Edition)
- ◆ Intel[®] Pentium[®] equivalent, processor speed: 1.0 GHz
- ◆ 512 Mbyte RAM
- ◆ 24 bit graphics card/monitor
- ◆ 1024 x 1280 pixel graphics (XVGA resolution monitor)

***syngo* Media Viewer Features & Benefits**

syngo Media Viewer is a comprehensive viewing application which can be recorded on a CD or DVD, enabling the distribution of clinical studies. With *syngo* Media Viewer, it is easy to share clinical images with referring physicians or export images for marketing and educational purposes. Simply select the patient's study and record the study to a disc with the Codonics Virtua Medical Disc Publisher, and the study is ready to be viewed on any PC.

The *syngo* Media Viewer allows you to visualize PET/CT or SPECT/CT images as well as stand-alone CT, MR, SPECT or PET studies. Fused studies are optimized for viewing with the images displayed in coronal, transaxial, and sagittal planes with a correlated MIP. Additionally, fused images are displayed in a format that allows blending between the PET or SPECT images and the CT. PET Standard Uptake Values (SUV) quantification is also included and can be calculated with this application.

With *syngo* Media Viewer, you have the ability to send entire studies to your referring physicians or medical facilities. Even users who do not have a fully DICOM-compliant image viewing workplace can view fused images using the tools provided by *syngo* Media Viewer. In addition to reading the report, your referring physicians can also navigate through a patient's entire scan.

***syngo* Media Viewer Features**

- ◆ View linked axial, coronal and sagittal views
- ◆ Navigate easily in three dimensions
- ◆ View MIP images correlated to coronal, transaxial, and sagittal views
- ◆ Blend fused images
- ◆ Quantify
 - ◆ SUV
 - ◆ Hounsfield Units
 - ◆ Voxel Intensity

Get it all with just one call
1-800-444-1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice.

Copyright © 2011 by Codonics, Inc. Printed in the U.S.A. Part No. 901-167-002 Rev. 03



17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com

Codonics Limited KK
New Shibaura Bldg. F1
1-3-11, Shibaura
Minato-ku, Tokyo, 105-0023 JAPAN
Phone: 81-3-5730-2297
Fax: 81-3-5730-2295