



Horizon® - Extremely rugged for any military application

Simplify Imaging

- **Print Versatility-** Horizon's unique blend of combining multiple sizes of diagnostic film, grayscale and color paper output in a single, compact device makes it the obvious choice for today's military and expandable ISO shelter environments. Offering superior image output, Horizon is the rugged, reliable, versatile solution for the military's harsh surroundings.
- **Low Cost Media-** Horizon's grayscale paper offers a low-cost solution for hardcopy documentation.

Rugged Design

- **Reliability-** Engineered to endure shock and vibration during military transportation, Horizon withstands the significant, prolonged exposure to rapid transitions in temperature with absolutely no degradation in image quality.
- **Simple Stowage-** Durable cabinet mounts in any ISO shelter or military environment while the docking tray secures the imager.

Improve Workflow

- **Easy to Site and Service -** Horizon's ultra-compact footprint allows it to be placed on a desk, counter or in any mobile application. It weighs less than 70 pounds and features standard power requirements with plug-and-play simplicity. Horizon's unique swap service and smart card technology eliminates the need for on-site service.



For more information contact Codonics at **800-444-1198** or **440-243-1198** and ask for our Military Sales Manager. Copyright © 2005 Codonics, Inc.