



Codonics Film

Codonics DirectVista® Clear Film cost-effectively delivers superior quality images in two sizes



14" x 17"
35 x 43 cm



8" x 10"
20 x 25 cm

"The quality is superior to the traditional film imagers we were previously using."

A. Joseph Borelli, M.D., MRI at Belfair - Bluffton, SC


- Completely digital, our unique direct thermal technology achieves image quality better or equal to "laser quality"
- Convenient dry printing eliminates the need for plumbing, special ventilation and siting requirements
- Multiple film sizes allow you to lower film costs by printing on the most economical size for your application

Codonics DirectVista® Thermal Grayscale Diagnostic Clear Film




Codonics patented imaging technology outputs superior diagnostic dry film with an unmatched image quality. DVC film is the ideal solution for the most demanding medical hardcopy applications. Completely eliminating space between lines by advancing the film in a smooth motion, the diagnostic output is precise and consistent. DVC film is available in 14" x 17" and 8" x 10", giving the user the option to scale the image to the appropriate film size, providing an alternative to print at a lower cost.


Depend on Codonics to deliver reliable, consistent images with the highest quality possible.



In thousands of side by side comparisons, board certified radiologists could not tell the difference between images printed on Codonics imagers from those printed on more expensive laser imagers. Codonics patented imaging technology outputs diagnostic film achieving unmatched versatility and image quality.



Dynamic Media Transport System (DMTS®) allows Codonics imagers to print while the media advances in a smooth motion eliminating space between lines.

<p>Specifications:</p> <ul style="list-style-type: none"> • Compatible with Horizon multi-media imagers • Continuous tone grayscale rendering • $D_{min} \leq 0.10$ OD, $D_{max} \geq 2.75$ OD (transmissive) • Thickness ~ 0.2mm 	<p>Available Sizes:</p> <p>Horizon Imagers 14" x 17" (35 x 43 cm) 8" x 10" (20 x 25 cm)</p> 
--	---



We bring the future into focus

Call Codonics today at **800.444.1198** or visit www.codonics.com for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: www.codonics.com/ip/patents Copyright©2006-2022 Codonics Inc.



17991 Englewood Drive
 Middleburg Heights, OH 44130 USA
 +1.440.243.1198
 +1.440.243.1334 Fax
 Email info@codonics.com
www.codonics.com