



Manufacturers Certificate of Quality

Hardcopy Thermal Reflective/Transparency Media

Product Code: 90 IWZ Radiographic Film **FDA Regulation:** 21CFR 892.1840
Japan KYOKA/Listing: 13BY6387

Types Covered: A-CVP/A4-CVP Horizon Media, 1600P-A/A4 CVP Media (Color Reflective)
 1600T-A/A4 CVT Media (Color Transparency Film)

CODONICS CVP/T	
Medical Classification/Medical Intended Use	
Gray Levels	256 level/16.7M Colors
Resolution	12.6 pixels/mm
Continuous Tone	YES
Abrasion Resistant	YES
Chemical/Water Resistant	YES
Stability/Consistency	Color Matching Controlled
Archiveability	>20 YEARS
Tear resistance	YES, Polypropylene Layer (P), PET(T)
Image Fade/Change	HEAT SENSITIVE/ NONE W/NORMAL VIEWING AND STORAGE
Repeatable Results	TESTED and PASS
Weight/Thickness	JIS P8143(weight) JIS 8118 (thickness)
Minimum Shelf Life	12 MO Under Suggested Storage Conditions
Density	> 2.2 MAXD
Linear Dynamic Range	<.25, > 2.2
STANDARDS	ANSI IT9.11/IT9.19
Harmonized Tariff Code:	3703.90.60 & 4811.90.8000
STORAGE AND TRANSPORT REQUIREMENTS	
Storage Unprinted Material	Extended Term Max 21 degrees C and 20-30% Relative Humidity
Nominal	Average Temp 18 C and 21 C
Extreme	Max 25 degrees C, Min 4 degrees C
Sustained	Max 28 degrees C for 3 successive weeks per year
Transport	Max 25 degrees C and Max 32 C for 1 week (total transport period)



**17991 Englewood Drive
 Middleburg Heights, Ohio 44130 USA
 Email: info@codonics.com
 www.codonics.com**