

### DirectVista Film

Codonics DirectVista Film is as a state-of-the-art direct-thermal film ideally suited for CT, MRI, ultrasound, functional imaging (nuclear medicine), CR, DR mammography, dental, portable X-ray, and other medical imaging applications.

### Specifications

- ◆ 320 DPI spatial resolution with Horizon Imagers
- ◆ 300 DPI spatial resolution with NP/EP Imagers
- ◆ Continuous tone grayscale rendering
- ◆ Blue Dmin < 0.25 OD, Dmax > 3.10 OD
- ◆ Clear Dmin < 0.10 OD, Dmax > 2.85 OD
- ◆ Thickness ~ 0.2mm

Size	Width (mm)	Length (mm)	Width (inches)	Length (inches)
14"x51"	354	1294	14.0	51.0
14"x36"	354	913	14.0	36.0
14"x17"	354	429	14.0	17.0
11"x14"	278	354	11.0	14.0
8"x10"	201	252	8.0	10.0

Size	Horizon (pixel image area)	NP/EP (pixel image area)
14"x51"	4322 x 15885	N/A
14"x36"	4322 x 11095	N/A
14"x17"	4322 x 5025	N/A
11"x14"	3376 x 4072	N/A
8"x10"	2406 x 2790	2280 x 2565

### Shelf Life

Greater than or equal to 26 months when stored at optimal storage conditions from date of manufacture (24 C° No relative humidity control).

### Archival


DirectVista Film complies with ISO 5466 standards for storage of photographic media. Accelerated testing shows that DirectVista Film is fully archival for a minimum of 20 years when stored at extended storage conditions (Max 21 C° and 20-30 % RH).

### Ordering Information

Imager	Media Size	Catalog Number
Horizon	14"x51"	1451-DVB
Horizon	14"x36"	1436-DVB
Horizon	14"x17"	1417-DVB, 1417-DVC
Horizon	11"x14"	1114-DVB, 1114-DVC
Horizon	8"x10"	810-DVB, 810-DVC
EP/NP-1660	8"x10"	1660T-810B, 1660T-810C

\* A suffix of "B" indicates blue base film and "C" indicates clear base film.

- ◆ DirectVista Film is shipped five boxes to the case (100 sheets/box for Horizon, 100 sheets/box for NP/EP).
- ◆ 14"x51" and 14"x36" DirectVista Film is shipped four boxes to the case (25 sheets/box).



Get it all with just one call  
**1-800-444-1198**

Copyright © 2005 by Codonics, Inc. All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Printed in the U.S.A. 901-130-001.c



Codonics Limited KK  
 New Shibaura Bldg. Fl. 1-3-11, Shibaura  
 Minato -ku, Tokyo, 105-0023 JAPAN.  
 81-3-5730-2297  
 81-3-5730-2295 Fax  
 Email info@codonics.com  
 www.codonics.com

17991 Englewood Drive  
 Middleburg Heights, OH 44130 USA  
 (440) 243-1198  
 (440) 243-1334 Fax  
 Email info@codonics.com  
 www.codonics.com