

Print.

Horizon® Multi-media Imager


Film, color and paper in one compact design obsoletes film-only imagers. The Horizon Multi-media Imager offers more versatility than any other imager, by simultaneously printing diagnostic quality film, color paper, and grayscale paper in one cost-effective, high performance device.

Horizon's multi-media alternatives allow you to show off the stunning output from any modality and market your advanced imaging equipment. Horizon is easy to use, requiring no operator intervention to go from film to color to grayscale printing. For CR/DR applications, our Horizon XL model prints digitally stitched images on one continuous film with "true-size" imaging up to 51" in length. Codonics film is daylight safe, easy to handle, and eliminates wet film processing.

Features/Benefits

- Horizon's all-in-one patented technology provides maximum versatility and cost savings
- Color paper shows off advanced post-processing applications / helps attract referral business
- Convenient, high quality, low-cost grayscale paper is room-light viewable, and perfect for patient files and referral copy with 20-year archivability
- Dry long film is ideal for scoliosis, long bone studies, and pediatric and adult spines
- Horizon offers the fastest time-to-first-print in its class, printing up to 100 films an hour and holding up to 300 sheets of film
- The Smart Card holds Horizon's configuration settings and is the basis for our revolutionary Sunrise Express "Swap" Warranty, a commitment that you're up and running within 24 hours of a service call
- Desktop size is perfect for tight spaces
- Meets industry's toughest shake and vibration standards. Ask about our custom mobile cabinets and mounting kits for mobile applications
- Exclusive custom job settings allow contrast and resolution to be uniquely set for each modality
- Validated with all major OEMs, PACS and workstations for seamless connectivity
- Models available to suit every need and budget

Codonics
edge-to-edge
archivable color
photographic print



VS






Dull paper
print prone to
fading and
damage



Horizon prints on several sizes of diagnostic film, grayscale paper and color paper -simultaneously, without operator intervention, all while sitting on a desktop.



14"x51", 14"x36"
Exclusive Dry Long
Media for Orthopaedics



14"x17", 11"x14", 8"x10", A, A4
Film, Grayscale and Color Paper