

# Horizon<sup>®</sup> Imager Software

**Release Notes**

---

## **Version 1.7.3 Release Notes**

Horizon version 1.7.3 software corrects two defects:

- ◆ “Edge blanking” has been disabled for 8x10 **DirectVista<sup>®</sup>** Film because of a rare problem where the outer edges of image data had the potential of being cropped. As a result, automatic centering has been disabled for 8x10 film. Image quality and printable area are unaffected. No other film sizes or media types are affected.
- ◆ An LPR protocol defect has been repaired allowing better compatibility with third-party framegrabber products.

In all other respects, version 1.7.3 is identical to version 1.7.2. See Version 1.7.2 Horizon Imager Software Release Notes (Codonics Part Number 901-084-003) for full details.

This software is compatible with all Horizon imagers. Any Horizon imager running software versions 1.4.x through 1.7.x can be upgraded to version 1.7.3.

---

*Get it all with just one call*  
**1-800-444-1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2000-2004 by Codonics, Inc. Printed in the U.S.A. Part No. 901-084-004 Rev.A.



17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
(440) 243-1198  
(440) 243-1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)

Codonics Limited KK  
New Shibaura Bldg. F1  
1-3-11, Shibaura  
Minato-ku, Tokyo, 105-0023 JAPAN  
Phone: 81-3-5730-2297  
Fax: 81-3-5730-2295