



**CODONICS**

*We bring the future into focus*

*Integrity<sup>™</sup> Medical Image  
Importer*

**IHE Integration Statement**

**Part Number: 900-409-001.03**

**Document / Software Version: 1.3.0**

**Date: 6/6/2008 4:20:00 PM**

<b>Revision</b>	<b>Date</b>	<b>Author</b>	<b>Description of Changes</b>
01	27 Mar 2008	Mike Lustig	Initial Release based on Virtua Integration statement by Glenn Burke
02	29 May 2008	Ross Goodman	Updated front page graphics
03	5 June 5, 2008	Ross Goodman	Added IRWF Profile

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Portable Data for Imaging	Portable Media Importer	N/A
IRWF compliant with Unscheduled Import and compliant with Scheduled Import WL entries. Not fully compliant for Scheduled Work Flow (SWF) at this time.	Portable Media Importer	N/A
<b>Internet address for vendor's IHE information:</b> <a href="http://www.codonics.com/ihe">www.codonics.com/ihe</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="http://www.codonics.com/Support/Integrity/pdf/Integrity_DICOM_Conformance_Statement.pdf">http://www.codonics.com/Support/Integrity/pdf/Integrity_DICOM_Conformance_Statement.pdf</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>